

## MULTIPOINT MAGNETIC STIRRER - HEXAMAG

'eltek' Multi Point Magnetic Stirrer christened as 'HexaMag' is a state of art magnetic stirrer which is excellent for stirring of light liquids. All the six stirring points can be independently controlled. The parameters of Speed, Timer, Heating Level and Direction can be set as per user requirement. HexaMag is extremely beneficial especially in laboratories & pharmaceutical industries which involve continuous stirring of bulk solutions.

### Technical Data

Model No.		
Dimension in mm (W X D X H)		
Weight		
Capacity	2 Liters on each stirring point	2 Liters on each stirring point
Stirring Points	6	6
Display		
Cabinet Moc	S.S. 304	S.S. 304
Max Speed	1200 RPM	1200 RPM
Timer		
Heater		
Rotation	Clockwise & Anticlockwise	Clockwise & Anticlockwise

**Multi Point Magnetic Stirrer with High Volumes are also offered.**

