

DEEP FREEZER (- 20°C)

MODEL NO.	EDF - 100	EDF - 200	EDF - 300	EDF - 400	EDF - 500	EDF - 600
Temperature	- 20° C	- 20° C	- 20° C	- 20° C	- 20° C	- 20° C

Salient Features and Technical Specifications for Deep Freezer are same as that of Quick Freezer except for working temperature. Also in some models environmentally friendly HFC refrigerant gas R-134 A is used which is CFC and HFC free

Optional Features:

- Manufacturing as per CE standards.
- GMP Models also available (Outer SS 304 1 mm thick & Inner SS 304 1 mm thick)
- Auto/ Manual Defrost system.
- Door Open Alarm
- Self Closing Doors
- Height adjustable shelves
- Battery Backup for display for 10 Hours along with Power Failure Alarm.
- 8/4 Channel temperature data scanner with sensors, interfacing RS 232 to RS 485 or USB Port
- 21 CFR Part 11 Compliant Software for Data Management.
- Mobile alert via GSM incase of system malfunction.
- IQ, OQ, PQ, DQ, calibration, validation and traceability documentation can be provided as per FDA protocols.

Calibration:

All models are calibrated with the help of master calibrator which is certified for its accuracy by Electronic Regional Test Laboratory Mumbai.

Quick Freezers & Deep Freezers can also be manufactured as per user-customised specification.

