

## Magnetic Stirrer without Hot Plate - AM-2D (ABS)



### Technical Features:

- The unit made of ABS Body with PU finish Coating fitted with M. S. base plate.
- The unit provided with LCD display for working parameters like Speed / Time
- Variable stirring capacity.
- Preset digital timer (1-99 Min.) or infinite mode.
- The unit fitted with PMDC motor for continuous operations.
- The unit provided with soft touch keypad for setting of Speed, Time
- Soft start and memory of last used program in case of power failure.
- Maximum speed – 1500 RPM
- Maximum Stirring Capacity -2 liters
- Stirring Paddle – 8 x 35 mm (Teflon Coated)

### Technical Data

| Type                    | AM2D        |
|-------------------------|-------------|
| Max. Speed RPM          | 1500        |
| Stirring Capacity Ltrs. | 2           |
| Dimension (WxDxH) mm    | 180x250x145 |
| Stirring Paddle ØxL mm  | 8x35        |
| Weight                  | 2.7kg.      |

## SMARTMIX SM 30A



'eltekt' **SMARTMIX** is the updated model of Vortex Mixer. The equipment is ideal for thorough quick mixing of liquids in Test tubes, small flasks and small bottles. The equipment provided with eccentric drive to a rubber buffer cup on the top. The container with the liquids to be mixed is held against the vibrating cup on top. The velocity of the agitation can be controlled by the speed of the vibrating cup and the pressure at which the container is held on to the vibrating cup. The unit is also provided with an adaptor with different micro tube holder to use it as a mini shaker. The adaptor can accommodate 6 micro tubes (Eppendorf tubes) of 0.4 ml, or 0.5 ml or 1.5 ml

### The Technical Specifications:

|                         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|-------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Housing                 | ABS moulded housing with poly Urethane coating mounted on steel base                                                                                                                                                                                                                                                                                                                                                                                                   |
| Stability mounting      | Vacuum rubber feet provided to keep the unit stable on table top to avoid skidding while in vibrating mode.                                                                                                                                                                                                                                                                                                                                                            |
| Drive Motor Assembly    | AC motor with eccentric drive mechanism to provide vibration to the rubber buffer cup                                                                                                                                                                                                                                                                                                                                                                                  |
| Speed Control           | Thyristorized Electronic speed control provided to facilitate stepless variable speed of vibration. Speed range 200 to 1600 RPM                                                                                                                                                                                                                                                                                                                                        |
| Auto / Manual mode      | <p>Provided with option of Auto / Manual Mode for operation.</p> <p><b>Auto Mode :</b> will start the vibration when the tube is pressed to the vibrating cup at the set speed and will stop when the tube is removed. The agitation varies depending on the set speed by the speed control and also depending on the pressure exerted on to the vibrating cup</p> <p><b>Manual Mode:</b> To manually switch on / off and will operate continuity till switch off.</p> |
| Adaptor for small tubes | The equipment is provided with a replaceable rubber adaptor that can accommodate micro tubes (Eppendorff tubes). 6 Nos of 0.4 ml, or 0.5 ml or 1.5 ml                                                                                                                                                                                                                                                                                                                  |
| Electric Power          | 220/230 Volts 50 Hz AC. The unit is provided with 5 Amps 3 pin plug and protective fuse.                                                                                                                                                                                                                                                                                                                                                                               |
| Dimensions-W x D x H    | 130 x 185 x 145 mm                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Weight                  | 3.3 kgs                                                                                                                                                                                                                                                                                                                                                                                                                                                                |

## elektrocrafft (India) pvt.ltd.

Plot no 35, Sector II, The Vasai Taluka Indl. Co-op. Estate Ltd., Gauripada, Vasai (E), Palghar - 401 208, Maharashtra India.

Tel: +91 - 250 - 2457033 / 7410041042 / 43 · Sales Direct : +91 9145693925

info@eltekindia.com · sales@eltekindia.com · www.eltekindia.com

Customercare : +91 8369469743 · customercare@eltekindia.com

· service@eltekindia.com

Branch Office : eltek.hyderabad@gmail.com

AUTHORISED DEALER

Dimensions and specifications given above are approximate and not binding. We reserve the right to make improvements modifications with out notice.