

BLOOD BAG TUBE SEALER

'eltek' Blood Bag Tube Sealer is made of finest materials which ensure efficient sealing of blood bag tubes at all times. Its advanced Radio Frequency Sealing System Ensures hermetic sealing minimizing hemolysis of blood. Tube segments can be easily separated after sealing. It is designed to protect itself from any mechanical or electrical hazards. Its ergonomic design makes it user friendly and flexible to be operated anywhere in blood banks or blood donation camps. It is light in weight which makes it easy to carry.



Technical Specifications :

MODEL NO.	ES-01	ES-02
Type	Bench Type	Bench Type with Battery Backup
Dimension in mm (W x D x H)	220 x 200 x 120	220 x 200 x 120
Unit Weight	4kg	5kg
Power Consumption	150VA	80VA
Power Supply	230VAC 50/60 Hz	230VAC 50/60 Hz
Output RF Frequency	40.65 MHz	40.65 MHz
Sealing Time (approx.)	1.5 sec.	1.5 sec.
Electrical Protection	5A Fuse	5A Fuse
Max. Tube Diameter	upto 6 mm	upto 6 mm
Tube Detection	Automatic	Automatic
Battery Back-Up	N.A.	Available
Capacity	500 Seals	800 Seals