


# BLOOD BANK WEIGHING SCALE

'eltek' Micro-controller based **Blood Bank Weighing Scale** is specially designed for weighing blood and blood components. It can measure accurate value of Blood components as per the specific gravity of the blood components.

## Salient Features :

- › Auto - conversion of weight to volume.
- › Different mode for Plasma Platelet and Whole Blood Weighing
- › LCD micro-processor based display.
- › Compact model.
- › Polycarbonate molded body.
- › Accuracy  $\pm 1\text{gm/ml}$ .



## Technical Specifications :

MODEL NO.	EWS-1D
Display	LCD
Accuracy	1 gm/1ml
Calibration	Auto
Transducer	Load Cell
Power Supply	12 VDC Adaptor
Feet	Adjustable
Weighing capacity	1 kg (1000gm)
Weight	1.5 kg
Dimension ( W x D x H )	195 x 250 x 100