

# BLOOD DONOR COUCH

'eltek' Micro Controller based **Blood Donor Couch** is designed for comfortable position for the Donor as per the hemodynamic principles. Donor's body position adjustable upright by using Remote control. Seat arrangement designed with consideration to maximize donor comfort. (The chair has special tray on which accessories like BP Apparatus, Blood collection Monitor can be placed conveniently) **Optional**.

## Salient Features :

- › Dual Heavy duty geared motors.
- › Castor wheel with locking ensure safety of the donor.
- › Micro control based controller.
- › Washable upholstery.
- › Main body duly powder coated.
- › Battery backup for uninterrupted function (Optional).
- › Interface facility for Blood collection monitor, BP Apparatus.
- › IV Stand (Optional)



## Technical Specifications :

MODEL	EDC-01	EDC-02
Power supply	230 VAC 50 Hz	230 VAC 50 Hz
Instrument Weight	60Kg	63Kg
Lifting Capacity	150Kg	150Kg
Power Consumption	100W	100W
No. of Motors	2	2
Length of seat	200cm	200cm
Width	68cm	68cm
Upholstery Thickness	50mm	50mm
Battery Backup	NA	2 to 3 Hrs