

## CLINICAL CENTRIFUGE

A range of centrifuges for general purpose use in Pathology Labs, Small Hospitals, Schools and Colleges. The centrifuge with 5 step speed regulator is provided with Fixed Angle Rotor Head. Optional Mechanical Timer of 60 minutes can be offered.

### Technical Data

Model	Unit	TC 815 A	TC 245 A
Max. Speed	RPM	4500	4000
Max. RCF	xg	2050	1650
Max. Tube Size	ml	15	5
Max. Capacity	ml	120	120
Width	mm	210	210
Depth	mm	235	235
Height	mm	230	255
Weight	kg	6	6
Connected Load	kVA	0.13	0.13



TC 815 A

- Vacutainer tubes of 5ml and 6-10ml can be accommodated with above Centrifuge