

Magnetic Stirrer without Hot Plate - AM-2D (ABS)



Technical Features:

- The unit made of ABS Body with PU finish Coating fitted with M. S. base plate.
- The unit provided with LCD display for working parameters like Speed / Time
- Variable stirring capacity.
- Preset digital timer (1-99 Min.) or infinite mode.
- The unit fitted with PMDC motor for continuous operations.
- The unit provided with soft touch keypad for setting of Speed, Time
- Soft start and memory of last used program in case of power failure.
- Maximum speed – 1500 RPM
- Maximum Stirring Capacity -2 liters
- Stirring Paddle – 8 x 35 mm (Teflon Coated)

Technical Data

Type	AM2D
Max. Speed RPM	1500
Stirring Capacity Ltrs.	2
Dimension (WxDxH) mm	180x250x145
Stirring Paddle ØxL mm	8x35
Weight	2.7kg.