

MEDICO CENTRIFUGE



TC 815 S



TC 415 ES

A range of centrifuges for general purpose use in Pathology Labs, Small Hospitals, Schools and Colleges. The centrifuge with 5 step speed regulator is provided with Horizontal (Swing Out) Rotor Head.

Technical Data

Model	Unit	TC 415 ES	TC 615 S	TC 815 S
Max. Speed	RPM	3700	3500	3500
Max. RCF	xg	1950	1900	1900
Max. Tube Size	ml	15	15	15
Max. Capacity	ml	60	90	120
Width	mm	295	310	310
Depth	mm	295	310	310
Height	mm	250	290	290
Weight	kg	6.00	6.50	6.50
Connected Load	kVA	0.11	0.13	0.13