

# MEDICOUCH

An ideal utility product for versatile use in Hospitals, Blood Banks, Testing Labs etc.

Suitable as Blood Donation couch, Dialysis Couch, Hospital Examination couch, Blood collection for lab testing - administration of injections as couch/ chair etc. The equipment with 3 motors and remote control for adjusting the positioning of the 3-part flexible couch to flat, reclining and sitting positions for the convenience of the user. The mounting lockable castor wheel ensure easy mobility as well as sturdy positioning as required.



## Technical Specification - eltek Medicouch EDC-03

01. Power Source	230 Volt A/C, 50Hz, 5 Amps
02 Lifting Capacity	150 Kg
03 Weight	90 Kg.
04 Initiation of motor movement	By remote control
05 No. of motor used	Three Motor
06 Head Rest	Adjustable
07 Manoeuvrability	Provided with lockable castor wheels
08 Dimensions	a) Leg Rest -660 mm
	b) Back Rest -900 mm
	c) Seat-540 mm
	d) Total Dimension- 2100 X 600 mm
	e) Width -600 mm
	f) Upholstery Thickness -100 mm
09 Color	Green With light grey
10 Manufacturing standards	ISO : 9001-2008
	ISO : 13485 -2016