

PLASMA EXTRACTOR / EXPRESSOR

'eltek' Plasma Extractor / Expressor is designed to extract blood components i.e. plasma from centrifuged blood bags. Its sturdy construction makes it durable and reliable, easy to use and portable. This unit is a manually or semi automatic operated instrument. The pressure is made of acrylic.



EPE-1



EPE-2

Technical Specifications of eltek Plasma Extractor :

MODEL NO.	EPE-1	EPI-2
Control System	Manual	Micro controller based controller
Pressure Plate	Transparent Acrylic	Transparent Acrylic
Material of construction	Mild steel powder coated	Mild steel powder coated
Size (mm)	L195x D250xH195	L195x D250xH220
Sensor	NA	IR
Alarm	NA	Audio visual
Clamping	NA	Motor activated
Power	NA	230 VAC 50Hz Single Phase