



PLATELET AGITATOR

'eltek' Platelet Agitator is designed to store platelet bags with good ventilation and continuous motion. Zink plated sliding trays can be removed for cleaning as well as easy removal and loading of bags in the tray. Heavy duty bearing for smooth operation. Interface facility with eltek platelet incubator for auto stop when door open.

Technical Specifications :

MODEL NO.	EPA-36	EPA-48
Capacity	36 bags	48 bags
Oscillation (Stroke)	70/stroke per minute	70/stroke per minute
Power Supply	230 VAC 50 Hz	230 VAC 50 Hz
Accuracy	± 5 stroke	± 5 stroke
Alarm (Optional)	Motion Failure	Motion Failure

