

VORTEX MIXER

The unit is meant for quick and thorough mixing of liquids in test tubes, small flasks or bottles. The mixer consists of an enamel finished housing containing the motor with an eccentric drive to a rubber buffer on the top. The tube is held against the vibrating buffer and the intensity of agitation controlled by varying the position in which the tube is held. The thyristorised speed controller enables adjustment of degree of oscillations. The unit is also provided with a separate adaptor for using it as a mini shaker for small tubes. The adaptor can accommodate 6 tubes of 0.4 ml, 6 tubes of 0.5 ml and 6 tubes of 1.5 ml, RPM 1500.

Technical Data

Type : VM 301	Connected Load : 30 Watts	RPM	Dimension (WxDxH) mm : 150x150x125	Weight : 3.1 kg
		1500		

The above unit operates on 220/230 V, 50 cycles, Ø A. C. Supply.

VORTEX MIXER

V 31 Adaptor for multiple tube holder - 6x0.4 ml, 6x0.5 ml, 6x1.5 ml.

V 32 Adaptor for tubes, bottles and flasks up to 250 ml.

